

Beckman Spinchron 15/R



Beckman GS-15/R Centrifuge

The Beckman GS-15 Series centrifuges are ideal for pelleting, binding studies, nucleic acid precipitation, immunoassays, protein extraction and concentration, and sample preparation. The GS-15 features a brushless induction drive motor which provides quiet, maintenance-free operation. Digital displays indicate actual speed, set speed, temperature and rotor imbalance with multiple acceleration and deceleration rates.



The Beckman GS-15 Series centrifuges use an automatic sensor blocks centrifuge operation until the door is locked. Once in operation, the sensor releases the door lock only when the centrifuge has stopped.

Specifications:

| | |
|-----------------------------------|--|
| Set Speed: | 0 to 15,300 rpm (in 100-rpm increments) or equivalent RCF |
| Speed Display: | Digital display indicates actual rotor speed +/- 1 rpm (actual RCF can be selected during operation) |
| Acceleration: | Three acceleration profiles |
| Deceleration: | Nine deceleration profiles |
| Ambient Temperature Range: | 10°C to 35°C |
| Temperature Setting: | -20°C to +40°C (1°C increments) |
| Time Display: | Time remaining in run (timed run) or infinity and elapsed time (continuous run) |
| Set Time: | To 9 hr 59 minutes or continuous |

Rotors and Capacities

| Rotor Description | Refrigerated | | Nonrefrigerated | | Max Capacity (mL) | Rotor Part Number |
|---------------------------------------|--------------|---------------|-----------------|---------------|-------------------|-------------------|
| | Max RPM | Max RCF (x g) | Max RPM | Max RCF (x g) | | |
| F2402H Fixed Angle | 15,300 | 21,460 | 14,500 | 19,280 | 24 x 1.0/1.5 | 361171 |
| F3602 Fixed Angle | 14,000 | 19,280 | 12,500 | 15,370 | 36 x 2.0/1.5 | 364600 |
| F0630 Fixed Angle | 15,300 | 20,410 | 14,500 | 18,340 | 6 x 30 | 361231 |
| F0850 Fixed Angle | 10,000 | 10,730 | 9,000 | 8,690 | 8 x 50 | 346640 |
| F0685 Fixed Angle | 8,000 | 7,010 | 8,000 | 7,010 | 6 x 85 | 364650 |
| F1010 Fixed Angle | 15,300 | 19,890 | 14,500 | 17,870 | 10 x 10 | 361221 |
| C0650 (Conical) Fixed Angle | 9,500 | 9,380 | 9,000 | 8,420 | 6 x 50 | 364670 |
| C1015 (Conical) Fixed Angle | 9,500 | 9,380 | 9,000 | 8,420 | 10 x 15 | 364680 |
| S4180 Swinging Buckets (4) | 5,500 | 5,450 | 4,500 | 3,650 | 4 x 180 | 361109 |
| S2096 Microliter | 3,000 | 1,100 | 3,000 | 1,100 | 96 x 1.5/1.8 | 361111 |