

## Silencer® S2450 Elite Series Centrifuge Parts List

Parts are sold individually. The quantity listed is the quantity used per unit.

Description	Quantity Needed	Part No.
Gasket, Lid	1 EA	SIL00452
Gasket, Motor	1 EA	SIL00456A
Hinge, Lid	2 EA	SIL00452
Latch, Lidlock	1 EA	SIL00643B
Lid	1 EA	SIL00660
Lidlock Assembly	1 EA	SIL00671B
Solenoid (part of SIL00671B)	1 EA	SIL00131
Switch, Lid Interlock (part of SIL00671B)	1 EA	SIL02456
Switch, Lock Interlock (part of SIL00671B)	1 EA	SIL02457
Motor Assembly	1 EA	OCM-18K
Mounting Bracket, Motor	2 EA	SIL00787
Nut, Rotor Securing	1 EA	SIL00791
PCB, Display/Control	1 EA	SIL00813D
PCB, Power	1 EA	35032
Relay, Drive / Brake	1 EA	SIL01103