



Centrifuges, Motors & Replacement Parts

Our motors meet or exceed the specifications of new motors.

The remanufacturing process includes all components necessary to bring the motor back to original specifications.

Our process includes new brushes, bearings, and any necessary motor housing repair or shaft replacement. A rewind armature with new commutator is used if needed.

All rebuilt armatures are balanced as well as dipped and baked to assure maximum performance. All motors are bench tested and most are instrument tested under a load.

All motors and spindles come with a 6-month full replacement warranty.

We sell motors and armatures 2 ways:

Exchanges are shipped the day ordered. Upon receiving the product we ask that you return your exchange in the same box within 15 days. All paper work and return shipping label are enclosed for your convenience. Your help returning Exchange Items on time is greatly appreciated!

Rebuild: You send your product to us for repair and receive the remanufactured product in return. Our turn-a-round is usually 1 day. Call for RGA# before sending any product to us! We recommend buying exchange motors when time is of utmost importance. It is not cost effective to request a rebuild over an exchange if you are going to ship priority both ways. With the exchange you can have your motor the next morning by 10:30AM.

