

Beckman Centrifuges

Spinchron 15 / Allegra 21



Display PCB Exchange
For 220V Non-Refrig. Units
OEM:329563
P/N:10001



Display PCB Exchange
For 120V Non-Refrig. Units
OEM:365745
P/N:10084



Power PCB Exchange
For 120V Units
OEM:365738
P/N:10129



Power PCB Exchange
For 220V Units
OEM:361357
P/N:10071



Lid Spring
Left 2.26"
OEM:363764
P/N:10055



Lid Spring
Right 2.26"
OEM:363765
P/N:10056



Tapered Sleeve
OEM:361328
P/N:10051



Motor Capseal
OEM:361330
P/N:10053



Touch Pad Exchange
OEM:365735
P/N:10068



Rotor ID Pickup
Fits the 120V Units
OEM:365752
P/N:10077



Tach Cable
Fits the 120V Units
OEM:365754
P/N:10078



Motor Exchange
Rebuild Price Available!
OEM:361348
P/N:OCM-021



115V Chassis Fan
OEM:961189
P/N:10090



220V Chassis Fan
OEM:961838
P/N:10103



120V Power Inlet
OEM:964643
P/N:10112



Tach PCB Exchange
Rectangular Style
OEM:361349
P/N:10049



Motor Mount Each
OEM:961803
P/N:10119



Rotor Tie Down
OEM:361367
P/N:10135



Motor Boot
OEM:364206
P/N:10064



New 120V Latch Exchange
OEM:365706
P/N:20122



120V Latch Exchange
OEM:365706
P/N:20123



220V Latch Exchange
OEM:961845
P/N:20127