

Sorvall

Cellwasher II



Motor Mount
OEM:90983
P/N:38080



5V Regulator
OEM:67078
P/N:38090



Capacitor
5600uF, 35V
OEM:67628
P/N:38095



Transformer
OEM:67635
P/N:38096



Flow Switch Assy.
OEM:04620
P/N:38024



Rubber Foot
OEM:68568
P/N:38099



O-Ring for Retainer
OEM:69722
P/N:38100



Fuse Kit
6.25A Slo-Blo (Qty. 5)
OEM:NA
P/N:38101



Saline Shield
OEM:04941
P/N:38116



Latch Solenoid
OEM:NA
P/N:38118



Cooling Fan
OEM:93006
P/N:38486



Lid Assy. Exchange
OEM:04660
P/N:38521



Pump Motor
OEM:04578
P/N:38020



New Bowl Assy. Exchange
OEM:12796
P/N:38049



Bowl Assy. Exchange
OEM:12796R
P/N:38002

Ozark Biomedical is a ISO9001:2008 Registered Company.

**Quality Parts and Services, along with Customer Satisfaction,
are Our Top Priorities.**